


Instructions for use

SonoCheck

SonoCheck is a dosimeter used for monitoring the Ultrasonic energy.



1. For the routine or validation test, SonoCheck is placed in the basket within the instrument load and the cycle to be tested is switched on. Regarding the number of SonoCheck to be used in different size bath, please refer to **the clean monitoring user guide**.
2. At the end of the cycle the vials are visually checked for the colour change from blue green to yellow. For a positive result SonoCheck need to change completely to yellow within the cycle time.
3. The time needed for the colour change can provide additional information regarding the ultrasonic energy.
4. The results are then recorded in the **SonoCheck log book**.
5. In case of unsatisfactory results please refer to the **troubleshooting guide**.

Note: To monitor the chemical cleaning efficacy of the process, the additional use of TOSI-LumCheck is recommended. Do not use SonoCheck for washer-disinfectors because of the high temperature during drying cycle.

STORAGE: Store SonoCheck-box at 2°C - 25°C and do not freeze
Keep away from light and heat

SHELF LIFE: See imprint

CONTENT: 30 tests

Marketing & Sales:

Manufacturer:
PEREG GMBH
Porschestr. 12
D-84478 Waldkraiburg
Germany